

兔抗 SLC6A3 多克隆抗体

中文名称：兔抗 SLC6A3 多克隆抗体

英文名称：Anti-SLC6A3 rabbit polyclonal antibody

别名：solute carrier family 6 member 3; DAT; DAT1; PKDYS; PKDYS1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SLC6A3

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a dopamine transporter which is a member of the sodium- and chloride-dependent neurotransmitter transporter family. The 3' UTR of this gene contains a 40 bp tandem repeat, referred to as a variable number tandem repeat or VNTR, which can be present in 3 to 11 copies. Variation in the number of repeats is associated with idiopathic epilepsy, attention-deficit hyperactivity disorder, dependence on alcohol and cocaine, susceptibility to Parkinson disease and protection against nicotine dependence.

Applications:	ELISA, IHC
Name of antibody:	SLC6A3
Immunogen:	Synthetic peptide of human SLC6A3
Full name:	solute carrier family 6 member 3
Synonyms:	DAT; DAT1; PKDYS; PKDYS1
SwissProt:	Q01959
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and Human brain
IHC Recommend dilution:	50-200



