

## 兔抗 ZBTB3 多克隆抗体

中文名称：兔抗 ZBTB3 多克隆抗体

英文名称：Anti-ZBTB3 rabbit polyclonal antibody

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ZBTB3

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	ZBTB3 (zinc finger and BTB domain containing 3) is a 574 amino acid protein that contains one BTB (POZ) domain and two C2H2-type zinc fingers. Localized to the nucleus, ZBTB3 is thought to play a role in transcriptional regulation events. The gene encoding ZBTB3 maps to human chromosome 11, which houses over 1,400 genes and comprises nearly 4% of the human genome.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ZBTB3
<b>Immunogen:</b>	Synthetic peptide of human ZBTB3
<b>Full name:</b>	zinc finger and BTB domain containing 3
<b>SwissProt:</b>	Q9H5J0

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and human colorectal cancer
<b>IHC Recommend dilution:</b>	20-100



