

## 兔抗 SLC6A6 多克隆抗体

中文名称：兔抗 SLC6A6 多克隆抗体

英文名称：Anti-SLC6A6 rabbit polyclonal antibody

别名：TAUT

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：SLC6A6

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

#### Background:

This gene encodes a multi-pass membrane protein that is a member of a family of sodium and chloride-ion dependent transporters. The encoded protein transports taurine and beta-alanine. Ther

	e is a pseudogene for this gene on chromosome 21. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	SLC6A6
<b>Immunogen:</b>	Synthesized peptide derived from internal of human SC6A6.
<b>Full name:</b>	solute carrier family 6 (neurotransmitter transporter), member 6
<b>Synonyms :</b>	TAUT
<b>SwissProt:</b>	P31641
<b>IHC positive control:</b>	Human lung carcinoma tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	70 kDa
<b>WB Positive control:</b>	COLO205 cells lysate
<b>WB Recommended dilution:</b>	500-3000

