

兔抗 SLC6A7 多克隆抗体

中文名称：兔抗 SLC6A7 多克隆抗体

英文名称： Anti-SLC6A7 rabbit polyclonal antibody

别名： PROT

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SLC6A7

反应种属： Human, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is a member of the gamma-aminobutyric acid (GABA) neurotransmitter gene family and encodes a high-affinity mammalian brain L-proline transporter protein. This transporter protein differs from other sodium-dependent plasma membrane carriers by its pharmacological specificity, kinetic properties, and ionic requirements.
Applications:	ELISA, IHC
Name of antibody:	SLC6A7
Immunogen:	Synthetic peptide of human SLC6A7
Full name:	solute carrier family 6 (neurotransmitter transporter),

	member 7
Synonyms :	PROT
SwissProt:	Q99884
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human brain and Human breast cancer
IHC Recommend dilution:	25-100

