

## 兔抗 KIAA1804 多克隆抗体

中文名称：兔抗 KIAA1804 多克隆抗体

英文名称：Anti-KIAA1804 rabbit polyclonal antibody

别名：MLK4; dJ862P8.3

抗原：KIAA1804

储存：冷冻（-20℃）避光

宿主：Rabbit

反应种属：Human

相关类别：一抗

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	Negative regulator of TLR4 signaling. Does not activate JN K1/MAPK8 pathway, p38/MAPK14, nor ERK2/MAPK1 pathways.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	KIAA1804
<b>Immunogen:</b>	Synthesized peptide derived from internal of human MAP3 KL4.
<b>Full name:</b>	mixed lineage kinase 4
<b>Synonyms :</b>	MLK4; dJ862P8.3
<b>SwissProt:</b>	Q5TCX8
<b>WB Predicted band size:</b>	114 kDa

<b>WB Positive control:</b>	HT-29 cells lysate
<b>WB Recommended dilution:</b>	500-3000

