

兔抗 SLC6A8 多克隆抗体

中文名称：兔抗 SLC6A8 多克隆抗体

英文名称： Anti-SLC6A8 rabbit polyclonal antibody

别名： CRT; CT1; CRTR; CTR5; CCDS1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SLC6A8

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a plasma membrane protein whose function is to transport creatine into and out of cells. Defects in this gene can result in X-linked creatine deficiency syndrome. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	SLC6A8
Immunogen:	Synthetic peptide of human SLC6A8
Full name:	solute carrier family 6 (neurotransmitter transporter), member 8

Synonyms :	CRT; CT1; CRTR; CTR5; CCDS1
SwissProt:	P48029
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and Human brain
IHC Recommend dilution:	25-100

