

## 兔抗 SLC6A9 多克隆抗体

中文名称：兔抗 SLC6A9 多克隆抗体

英文名称：Anti-SLC6A9 rabbit polyclonal antibody

别名：solute carrier family 6 member 9; GLYT1; GCENSG

抗原：SLC6A9

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The amino acid glycine acts as an inhibitory neurotransmitter in the central nervous system. The protein encoded by this gene is one of two transporters that stop glycine signaling by removing it from the synaptic cleft.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SLC6A9
<b>Immunogen:</b>	Synthetic peptide of human SLC6A9
<b>Full name:</b>	solute carrier family 6 member 9
<b>Synonyms:</b>	GLYT1; GCENSG
<b>SwissProt:</b>	P48067
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human colorectal cancer and Human cervical cancer
<b>IHC Recommend dilution:</b>	30-150



