

## 兔抗 SLC7A1 多克隆抗体

中文名称：兔抗 SLC7A1 多克隆抗体

英文名称：Anti-SLC7A1 rabbit polyclonal antibody

别名：solute carrier family 7 member 1; ERR; ATRC1; CAT-1; HCAT1; REC1L

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SLC7A1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	High-affinity, low capacity permease involved in the transport of the cationic amino acids (arginine, lysine and ornithine) in non-hepatic tissues. May also function as an ectopic retroviral leukemia receptor.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SLC7A1
<b>Immunogen:</b>	Synthetic peptide of human SLC7A1
<b>Full name:</b>	solute carrier family 7 member 1
<b>Synonyms:</b>	ERR; ATRC1; CAT-1; HCAT1; REC1L
<b>SwissProt:</b>	P30825
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-100

