

## 兔抗 PLPP1 多克隆抗体

中文名称：兔抗 PLPP1 多克隆抗体

英文名称： Anti-PLPP1 rabbit polyclonal antibody

别名： phospholipid phosphatase 1; LPP1; PAP2; LLP1a; PAP

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PLPP1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a member of the phosphatidic acid phosphatase (PAP) family. PAPs convert phosphatidic acid to diacylglycerol, and function in synthesis of glycerolipids and in phospholipase D-mediated signal transduction. This enzyme is an integral membrane glycoprotein that plays a role in the hydrolysis and uptake of lipids from extracellular space. Alternate splicing results in multiple transcript variants of this gene. [provided by RefSeq, May 2013]
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PLPP1

<b>Immunogen:</b>	Synthetic peptide of human PLPP1
<b>Full name:</b>	phospholipid phosphatase 1
<b>Synonyms:</b>	LPP1; PAP2; LLP1a; PAP-2a; PPAP2A
<b>SwissProt:</b>	O14494
<b>IHC positive control:</b>	Human liver cancer and Human lung cancer
<b>IHC Recommend dilution:</b>	50-200

