

兔抗 ZC2HC1A 多克隆抗体

中文名称：兔抗 ZC2HC1A 多克隆抗体

英文名称：Anti-ZC2HC1A rabbit polyclonal antibody

别名：zinc finger C2HC-type containing 1A; CGI-62; C8orf70; FAM164A

相关类别：一抗

抗原：ZC2HC1A

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, IHC
Name of antibody:	ZC2HC1A
Immunogen:	Synthetic peptide of human ZC2HC1A
Full name:	zinc finger C2HC-type containing 1A
Synonyms:	CGI-62; C8orf70; FAM164A
SwissProt:	Q96GY0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human prostate cancer
IHC Recommend dilution:	20-100



