

## 兔抗 PLXNA1 多克隆抗体

中文名称: 兔抗 PLXNA1 多克隆抗体

英文名称: Anti-PLXNA1 rabbit polyclonal antibody

别名: NOV; NOVP; PLXN1; PLEXIN-A1

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: PLXNA1

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Plexins are a family of large, transmembrane receptors for multiple classes of semaphorins in vertebrates. Plexins are widely expressed, and regions of their extracellular domain are homologous to both scatter factor receptors and semaphorin domains.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PLXNA1
<b>Immunogen:</b>	Synthetic peptide of human PLXNA1
<b>Full name:</b>	plexin A1
<b>Synonyms :</b>	NOV; NOVP; PLXN1; PLEXIN-A1

<b>SwissProt:</b>	Q9UIW2
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human prostate cancer and Human brain
<b>IHC Recommend dilution:</b>	50-200

