

兔抗 PMS1 多克隆抗体

中文名称: 兔抗 PMS1 多克隆抗体

英文名称: Anti-PMS1 rabbit polyclonal antibody

别名: DKFZp781M0253; FLJ98259; HNPCC3; PMSL1; hPMS1; PMS1

相关类别: 一抗

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

抗原: PMS1

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a protein belonging to the DNA mismatch repair mutL/hexB family. This protein is thought to be involved in the repair of DNA mismatches, and it can form heterodimers with MLH1, a known DNA mismatch repair protein. Mutations in this gene cause hereditary nonpolyposis colorectal cancer type 3 (HNPCC3) either alone or in combination with mutations in other genes involved in the HNPCC phenotype, which is also known as Lynch syndrome.
Applications:	WB
Name of antibody:	PMS1

Immunogen:	Fusion protein of human PMS1
Full name:	PMS1 homolog 1, mismatch repair system component
Synonyms:	MLH2; PMSL1; hPMS1; HNPCC3
SwissProt:	P54277
WB Predicted band size:	106 kDa
WB Positive control:	HeLa, Jurkat, Molt-4 and HepG2 cells
WB Recommended dilution:	500-2000

