

兔抗 SLC9A7 多克隆抗体

中文名称：兔抗 SLC9A7 多克隆抗体

英文名称：Anti-SLC9A7 rabbit polyclonal antibody

别名：NHE7; NHE-7; SLC9A6

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SLC9A7

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a sodium and potassium/ proton antiporter that is a member of the solute carrier family 9 protein family. The encoded protein is primarily localized to the trans-Golgi network and is involved in maintaining pH homeostasis in organelles along the secretory and endocytic pathways. This protein may enhance cell growth of certain breast tumors. This gene is part of a gene cluster on chromosome Xp11.23. A pseudogene of this gene is found on chromosome 12. Alternate splicing results in multi

	ple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	SLC9A7
Immunogen:	Synthetic peptide of human SLC9A7
Full name:	solute carrier family 9, subfamily A (NHE7, cation/proton antiporter 7), member 7
Synonyms :	NHE7; NHE-7; SLC9A6
SwissProt:	Q96T83
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	50-200

