

兔抗 SLC9A9 多克隆抗体

中文名称：兔抗 SLC9A9 多克隆抗体

英文名称： Anti-SLC9A9 rabbit polyclonal antibody

别名： NHE9; AUTS16

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SLC9A9

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a sodium/proton exchanger that is a member of the solute carrier 9 protein family. The encoded protein localizes to the late recycling endosomes and may play an important role in maintaining cation homeostasis. Mutations in this gene are associated with autism susceptibility 16 and attention-deficit/hyperactivity disorder.
Applications:	ELISA, IHC
Name of antibody:	SLC9A9
Immunogen:	Synthetic peptide of human SLC9A9
Full name:	solute carrier family 9, subfamily A (NHE9, cation pro

	ton antiporter 9), member 9
Synonyms :	NHE9; AUTS16
SwissProt:	Q8IVB4
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human breast cancer
IHC Recommend dilution:	25-100

