

兔抗 PNLIPRP1 多克隆抗体

中文名称：兔抗 PNLIPRP1 多克隆抗体

英文名称：Anti-PNLIPRP1 rabbit polyclonal antibody

别名：PLRP1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PNLIPRP1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

WB Recommended dilution:	200-1000
WB Positive control:	K562 and 293T cell lysates
WB Predicted band size:	52 kDa
IHC Recommend dilution:	10-50
IHC positive control:	Human brain
ELISA Recommended dilution:	5000-10000
SwissProt:	P54315
Synonyms:	PLRP1

Full name:	pancreatic lipase-related protein 1
Immunogen:	Fusion protein of human PNLIPRP1
Name of antibody:	PNLIPRP1
Applications:	ELISA, WB, IHC
Background:	May function as inhibitor of dietary triglyceride digestion. Lacks detectable lipase activity towards triglycerides, diglycerides, phosphatidylcholine, galactolipids or cholesterol esters (in vitro) (By similarity).



