

兔抗 ZC3H13 多克隆抗体

中文名称：兔抗 ZC3H13 多克隆抗体

英文名称：Anti-ZC3H13 rabbit polyclonal antibody

别名：KIAA0853

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：ZC3H13

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	2 isoforms of the human protein are produced by alternative splicing.
Applications:	WB
Name of antibody:	ZC3H13
Immunogen:	Synthesized peptide derived from internal of human ZC3H13.
Full name:	zinc finger CCCH-type containing 1
Synonyms :	KIAA0853
SwissProt:	Q5T200
WB Predicted band size:	197 kDa
WB Positive control:	HepG2 cells, A549 cells and Jurkat cells lysates

WB Recommended dilution: 500-3000

