

## 兔抗 SLC11A2 多克隆抗体

中文名称：兔抗 SLC11A2 多克隆抗体

英文名称：Anti-SLC11A2 rabbit polyclonal antibody

别名：DCT1; DMT1; NRAMP2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SLC11A2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the solute carrier family 11 protein family. The product of this gene transports divalent metals and is involved in iron absorption. Mutations in this gene are associated with hy
--------------------	--

	<p>pochromic microcytic anemia with iron overload. A related solute carrier family 11 protein gene is located on chromosome 2 . Multiple transcript variants encoding different isoforms have been found for this gene.</p>
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SLC11A2
<b>Immunogen:</b>	Synthetic peptide of human SLC11A2
<b>Full name:</b>	solute carrier family 11 (proton-coupled divalent metal ion transporter), member 2
<b>Synonyms :</b>	DCT1; DMT1; NRAMP2
<b>SwissProt:</b>	P49281
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	50-200

