

兔抗 PNPLA6 多克隆抗体

中文名称: 兔抗 PNPLA6 多克隆抗体

英文名称: Anti-PNPLA6 rabbit polyclonal antibody

别名: NTE; sws; BNHS; SPG39; NTEMND; iPLA2delta

抗原: PNPLA6

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:

This gene encodes a phospholipase that deacetylates intracellular phosphatidylcholine to produce glycerophosphocholine. It is thought to function in neurite outgrowth and process elongation during neuronal differentiation. The protein is anchored to the cytoplasmic face of the endoplasmic reticulum in both neurons and non-neuronal cells. Mutations in this gene result in autosomal recessive spastic paraplegia, and the protein is the target for neurodegeneration induced by organophosphorus compounds and chemical warfare agents. Multiple transcript variants encoding differen

| | |
|---------------------------------|--|
| | t isoforms have been found for this gene. |
| Applications: | WB |
| Name of antibody: | PNPLA6 |
| Immunogen: | Synthesized peptide derived from internal of human PNPLA6. |
| Full name: | patatin-like phospholipase domain containing 6 |
| Synonyms : | NTE; sws; BNHS; SPG39; NTEMND; iPLA2delta |
| SwissProt: | Q8IY17 |
| WB Predicted band size: | 150 kDa |
| WB Positive control: | COLO205 cells lysate |
| WB Recommended dilution: | 500-3000 |

