

兔抗 ZDHHC2 多克隆抗体

中文名称：兔抗 ZDHHC2 多克隆抗体

英文名称：Anti-ZDHHC2 rabbit polyclonal antibody

别名：DHHC2; ZNF372

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：ZDHHC2

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Palmitoyltransferase specific for GAP43 and DLG4/PS D95 By similarity.
Applications:	WB
Name of antibody:	ZDHHC2
Immunogen:	Synthesized peptide derived from internal of human ZDHHC2.
Full name:	zinc finger, DHHC-type containing 2
Synonyms :	DHHC2; ZNF372
SwissProt:	Q9UIJ5
WB Predicted band size:	42 kDa
WB Positive control:	COS7 cells lysate

WB Recommended dilution:

500-3000

