

兔抗 SLC12A5 多克隆抗体

中文名称：兔抗 SLC12A5 多克隆抗体

英文名称： Anti-SLC12A5 rabbit polyclonal antibody

别名： solute carrier family 12 member 5; KCC2; EIG14; hKCC2; EIEE34

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SLC12A5

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	K-Cl cotransporters are proteins that lower intracellular chloride concentrations below the electrochemical equilibrium potential. The protein encoded by this gene is an integral membrane K-Cl cotransporter that can function in either a net efflux or influx pathway, depending on the chemical concentration gradients of potassium and chloride. The encoded protein can act as a homomultimer, or as a heteromultimer with other K-Cl cotransporters, to maintain chloride homeostasis in neurons. Alternative splicing results in two transcript variants encoding different isoforms.
Applications:	ELISA, IHC

Name of antibody:	SLC12A5
Immunogen:	Synthetic peptide of human SLC12A5
Full name:	solute carrier family 12 member 5
Synonyms:	KCC2; EIG14; hKCC2; EIEE34
SwissProt:	Q9H2X9
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	30-150

