

## 兔抗 POC1A 多克隆抗体

中文名称：兔抗 POC1A 多克隆抗体

英文名称： Anti-POC1A rabbit polyclonal antibody

别名： POC1 centriolar protein A; PIX2; SOFT; WDR51A

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： POC1A

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	POC1 proteins contain an N-terminal WD40 domain and a C-terminal coiled coil domain and are part of centrosomes. They play an important role in basal body and cilia formation. This gene encodes one of the two POC1 proteins found in humans. Mutations in this gene result in short stature, onychodysplasia, facial dysmorphism, and hypotrichosis (SOFT) syndrome.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	POC1A
<b>Immunogen:</b>	Fusion protein of human POC1A

<b>Full name:</b>	POC1 centriolar protein A
<b>Synonyms:</b>	PIX2; SOFT; WDR51A
<b>SwissProt:</b>	Q8NBT0
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	25-50

