

## 兔抗 SLC13A3 多克隆抗体

中文名称：兔抗 SLC13A3 多克隆抗体

英文名称： Anti-SLC13A3 rabbit polyclonal antibody

别名： NADC3; SDCT2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SLC13A3

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

#### Background:

Mammalian sodium-dicarboxylate cotransporters transport succinate and other Krebs cycle intermediates. They fall into 2 categories based on their substrate affinity: low affinity and high affinity. Both the low- and high-affinity transporters play an important role in the handling of citrate by the kidneys. The protein encoded by this gene represents the high-affinity form. Alternatively spliced transcript variants encoding different isoforms have been found for this gene, although the full-length nature of some of them have not been characterized yet.

<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	SLC13A3
<b>Immunogen:</b>	Synthetic peptide of human SLC13A3
<b>Full name:</b>	solute carrier family 13 (sodium-dependent dicarboxylate transporter), member 3
<b>Synonyms :</b>	NADC3; SDCT2
<b>SwissProt:</b>	Q8WWT9
<b>ELISA Recommended dilution:</b>	1000-2000
<b>WB Predicted band size:</b>	67 kDa
<b>WB Positive control:</b>	Human plasma and placenta tissue
<b>WB Recommended dilution:</b>	200-1000

