

兔抗 PODXL 多克隆抗体

中文名称：兔抗 PODXL 多克隆抗体

英文名称：Anti-PODXL rabbit polyclonal antibody

别名：PC, PCLP, Gp200, PCLP-1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PODXL

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a member of the sialomucin protein family. The encoded protein was originally identified as an important component of glomerular podocytes. Podocytes are highly differentiated epithelial cells with interdigitating foot processes covering the outer aspect of the glomerular basement membrane. Other biological activities of the encoded protein include: binding in a membrane protein complex with Na⁺/H⁺ exchanger regulatory factor to intracellular cytoskeletal elements, playing a role in hematopoietic cell differentiation, and being expressed

	in vascular endothelium cells and binding to L-selectin.
Applications:	ELISA, IHC
Name of antibody:	PODXL
Immunogen:	Synthetic peptide of human PODXL
Full name:	Podocalyxin-like
Synonyms :	PC, PCLP, Gp200, PCLP-1
SwissProt:	O00592
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human colon cancer and Human ovarian cancer
IHC Recommend dilution:	25-100

