

兔抗 SLC15A2 多克隆抗体

中文名称：兔抗 SLC15A2 多克隆抗体

英文名称： Anti-SLC15A2 rabbit polyclonal antibody

别名： PEPT2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SLC15A2

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The mammalian kidney expresses a proton-coupled peptide transporter that is responsible for the absorption of small peptides, as well as beta-lactam antibiotics and other peptide-like drugs, from the tubular filtrate. This transporter, SLC15A2, belongs to the same gene family as SLC15A1 (MIM 600544), the proton-coupled peptide transporter found in the small intestine
Applications:	ELISA, IHC
Name of antibody:	SLC15A2
Immunogen:	Synthetic peptide of human SLC15A2
Full name:	solute carrier family 15 (H ⁺ /peptide transporter), member

	2
Synonyms :	PEPT2
SwissProt:	Q16348
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human lung cancer
IHC Recommend dilution:	50-200

