

## 兔抗 SLC16A1 多克隆抗体

中文名称：兔抗 SLC16A1 多克隆抗体

英文名称： Anti-SLC16A1 rabbit polyclonal antibody

别名： MCT, HHF7, MCT1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SLC16A1

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a proton-linked monocarboxylate transporter that catalyzes the movement of many monocarboxylates, such as lactate and pyruvate, across the plasma membrane. Mutations in this gene are associated with erythrocyte lactate transporter defect. Alternatively spliced transcript variants have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SLC16A1
<b>Immunogen:</b>	Synthetic peptide of human SLC16A1

<b>Full name:</b>	Solute carrier family 16, member 1 (monocarboxylic acid transporter 1)
<b>Synonyms :</b>	MCT, HHF7, MCT1
<b>SwissProt:</b>	P53985
<b>ELISA Recommended dilution:</b>	2000-10000
<b>IHC positive control:</b>	Human liver cancer and Human ovarian cancer
<b>IHC Recommend dilution:</b>	50-200

