

兔抗 POLB 多克隆抗体

中文名称: 兔抗 POLB 多克隆抗体

英文名称: Anti-POLB rabbit polyclonal antibody

别名: DNA polymerase beta

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: POLB

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a DNA polymerase involved in base excision and repair, also called gap-filling DNA synthesis. The encoded protein, acting as a monomer, is normally found in the cytoplasm, but it translocates to the nucleus upon DNA damage. Several transcript variants of this gene exist, but the full-length nature of only one has been described to date. [provided by RefSeq, Sep 2011]
Applications:	ELISA, WB, IHC
Name of antibody:	POLB

Immunogen:	Fusion protein of human POLB
Full name:	DNA polymerase beta
SwissProt:	P06746
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human gastric cancer
IHC Recommend dilution:	100-300
WB Predicted band size:	38 kDa
WB Positive control:	Jurkat cell lysate
WB Recommended dilution:	500-2000



