

兔抗 ZDHHC15 多克隆抗体

中文名称: 兔抗 ZDHHC15 多克隆抗体

英文名称: Anti-ZDHHC15 rabbit polyclonal antibody

别名: MRX91

储存: 冷冻 (-20℃) 避光

抗原: ZDHHC15

宿主: Rabbit

反应种属: Human Mouse Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	The protein encoded by this gene belongs to the DHHC palmitoyltransferase family. Mutations in this gene are associated with mental retardatio X-linked type 91 (MRX91). Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	WB
Name of antibody:	ZDHHC15
Immunogen:	Synthesized peptide derived from C-terminal of human ZDHHC15.
Full name:	zinc finger, DHHC-type containing 15
Synonyms :	MRX91

SwissProt:	Q96MV8
WB Predicted band size:	39 kDa
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000

