

兔抗 SLC16A6 多克隆抗体

中文名称：兔抗 SLC16A6 多克隆抗体

英文名称： Anti-SLC16A6 rabbit polyclonal antibody

别名： MCT6; MCT7

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SLC16A6

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Proton-linked monocarboxylate transporter. Catalyzes the rapid transport across the plasma membrane of many monocarboxylates such as lactate, pyruvate, branched-chain oxo acids derived from leucine, valine and isoleucine, and the ketone bodies acetoacetate, beta-hydroxybutyrate and acetate.
Applications:	ELISA, IHC
Name of antibody:	SLC16A6
Immunogen:	Synthetic peptide of human SLC16A6
Full name:	solute carrier family 16 member 6
Synonyms:	MCT6; MCT7

SwissProt:	O15403
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer
IHC Recommend dilution:	50-100

