

## 兔抗 KIT(Phospho-Tyr703) 多克隆抗体

中文名称:	兔抗 KIT(Phospho-Tyr703) 多克隆抗体
英文名称:	Anti-KIT(Phospho-Tyr703) rabbit polyclonal antibody
别名:	PBT; SCFR; C-Kit; CD117
相关类别:	一抗
储存:	冷冻 (-20℃) 避光
宿主:	Rabbit
抗原:	KIT(Phospho-Tyr703)
反应种属:	Human
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes the human homolog of the proto-oncogene c-kit. C-kit was first identified as the cellular homolog of the feline sarcoma viral oncogene v-kit. This protein is a type 3 transmembrane receptor for MGF (mast cell growth factor, also known as stem cell factor). Mutations in this gene are associated with gastrointestinal stromal tumors, mast cell disease, acute myelogenous leukemia, and piebaldism. Multiple transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	KIT(Phospho-Tyr703)

<b>Immunogen:</b>	Peptide sequence around phosphorylation site of Tyrosine 703 (A-L-Y(p)-K-N) derived from Human KIT.
<b>Full name:</b>	v-kit Hardy-Zuckerman 4 feline sarcoma viral oncogene homolog
<b>Synonyms :</b>	PBT; SCFR; C-Kit; CD117
<b>SwissProt:</b>	P10721
<b>WB Predicted band size:</b>	110 kDa
<b>WB Positive control:</b>	HepG2 cells lysates treated with EGF
<b>WB Recommended dilution:</b>	500-1000

