

兔抗 POLL 多克隆抗体

中文名称：兔抗 POLL 多克隆抗体

英文名称： Anti-POLL rabbit polyclonal antibody

别名： BETAN; POLKAPPA

抗原： POLL

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a DNA polymerase. DNA polymerases catalyze DNA-template-directed extension of the 3'-end of a DNA strand. This particular polymerase, which is a member of the X family of DNA polymerases, likely plays a role in non-homologous end joining and other DNA repair processes. Alternatively spliced transcript variants have been described.
Applications:	WB
Name of antibody:	POLL
Immunogen:	Synthesized peptide derived from internal of human DNA P

	olymerase 位.
Full name:	polymerase (DNA directed), lambda
Synonyms :	BETAN; POLKAPPA
SwissProt:	Q9UGP5
WB Predicted band size:	63 kDa
WB Positive control:	293 cells lysate
WB Recommended dilution:	500-3000

