

## 兔抗 POLR1B 多克隆抗体

中文名称：兔抗 POLR1B 多克隆抗体

英文名称： Anti-POLR1B rabbit polyclonal antibody

别名： RPA2; RPA135; Rpo1-2

抗原： POLR1B

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	Eukaryotic RNA polymerase I (pol I) is responsible for the transcription of ribosomal RNA (rRNA) genes and production of rRNA, the primary component of ribosomes. Pol I is a multisubunit enzyme composed of 6 to 14 polypeptides, depending on the species. Most of the mass of the pol I complex derives from the 2 largest subunits, Rpa1 and Rpa2 in yeast.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	POLR1B
<b>Immunogen:</b>	Synthesized peptide derived from C-terminal of human POL

	R1B.
<b>Full name:</b>	polymerase (RNA) I polypeptide B, 128kDa
<b>Synonyms :</b>	RPA2; RPA135; Rpo1-2
<b>SwissProt:</b>	Q9H9Y6
<b>WB Predicted band size:</b>	128 kDa
<b>WB Positive control:</b>	Jurkat and NIH/3T3 cells lysates
<b>WB Recommended dilution:</b>	500-3000

