

兔抗 KLF6 多克隆抗体

中文名称：兔抗 KLF6 多克隆抗体

英文名称：Anti-KLF6 rabbit polyclonal antibody

别名：GBF; ZF9; BCD1; CBA1; CPBP; PAC1; ST12; COPEB

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：KLF6

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the Kruppel-like family of transcription factors. The zinc finger protein is a transcriptional activator, and functions as a tumor suppressor. Multiple transcript variants encoding different isoforms have been found for this gene, some of which are implicated in carcinogenesis.
Applications:	ELISA, IHC
Name of antibody:	KLF6
Immunogen:	Synthetic peptide of human KLF6
Full name:	Kruppel-like factor 6
Synonyms:	GBF; ZF9; BCD1; CBA1; CPBP; PAC1; ST12; COPEB

SwissProt:	Q99612
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

