

## 兔抗 POLR2D 多克隆抗体

中文名称: 兔抗 POLR2D 多克隆抗体

英文名称: Anti-POLR2D rabbit polyclonal antibody

别名: RBP4; RPB16; HSRBP4; HSRPB4

相关类别: 一抗

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

抗原: POLR2D

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes the fourth largest subunit of RNA polymerase II, the polymerase responsible for synthesizing messenger RNA in eukaryotes. In yeast, this polymerase subunit is associated with the polymerase under suboptimal growth conditions and may have a stress protective role. A sequence for a ribosomal pseudogene is contained within the 3' untranslated region of the transcript from this gene. [provided by RefSeq, Jul 2008]
<b>Applications:</b>	WB
<b>Name of antibody:</b>	POLR2D

<b>Immunogen:</b>	Fusion protein of human POLR2D
<b>Full name:</b>	polymerase (RNA) II (DNA directed) polypeptide D
<b>Synonyms :</b>	RBP4; RPB16; HSRBP4; HSRPB4
<b>SwissProt:</b>	O15514
<b>WB Predicted band size:</b>	16 kDa
<b>WB Positive control:</b>	Heart tissue
<b>WB Recommended dilution:</b>	500-2000

