

兔抗 KLF15 多克隆抗体

中文名称：兔抗 KLF15 多克隆抗体

英文名称：Anti-KLF15 rabbit polyclonal antibody

别名：KKLF

抗原：KLF15

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Krüppel-like factor 15 is a protein that in humans is encoded by the KLF15 gene in the Krüppel-like factor family. KLF15, a kidney-enriched Krüppel-like factor, is a transcriptional activator that binds the CLCNKA promoter. KLF15, which is a nuclear protein, is expressed primarily in liver, heart, skeletal muscle and kidney tissues but is not detected in lymphoid tissues or bone marrow. It is an important regulator of GLUT4 in both adipose and muscle tissues.
Applications:	ELISA, IHC
Name of antibody:	KLF15
Immunogen:	Fusion protein of human KLF15

Full name:	Kruppel-like factor 15
Synonyms :	KKLF
SwissProt:	Q9UIH9
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and human breast cancer
IHC Recommend dilution:	25-100

