

兔抗 POLR2I 多克隆抗体

中文名称：兔抗 POLR2I 多克隆抗体

英文名称：Anti-POLR2I rabbit polyclonal antibody

别名：RPB9; hRPB14.5

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：POLR2I

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a subunit of RNA polymerase II, the polymerase responsible for synthesizing messenger RNA in eukaryotes. This subunit, in combination with two other polymerase subunits, forms the DNA binding domain of the polymerase, a groove in which the DNA template is transcribed into RNA. The product of this gene has two zinc finger motifs with conserved cysteines and the subunit does possess zinc binding activity.
Applications:	WB
Name of antibody:	POLR2I
Immunogen:	Fusion protein of human POLR2I

Full name:	polymerase (RNA) II (DNA directed) polypeptide I
Synonyms :	RPB9; hRPB14.5
SwissProt:	P36954
WB Predicted band size:	15 kDa
WB Positive control:	LO2 cells
WB Recommended dilution:	500-2000

