

兔抗 KLHDC2 多克隆抗体

中文名称：兔抗 KLHDC2 多克隆抗体

英文名称：Anti-KLHDC2 rabbit polyclonal antibody

别名：kelch domain containing 2; LCP; HCLP1; HCLP-1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：KLHDC2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Represses CREB3-mediated transcription by interfering with CREB3-DNA binding.
Applications:	ELISA, IHC
Name of antibody:	KLHDC2
Immunogen:	Fusion protein of human KLHDC2
Full name:	kelch domain containing 2
Synonyms:	LCP; HCLP1; HCLP-1
SwissProt:	Q9Y2U9
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer

IHC Recommend dilution: 30-150

