

兔抗 POLR3B 多克隆抗体

中文名称：兔抗 POLR3B 多克隆抗体

英文名称： Anti-POLR3B rabbit polyclonal antibody

别名： C128; HLD8; RPC2; INMAP

抗原： POLR3B

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属 Human Mouse

相关类别 一抗

标记物 Unconjugate

克隆类型 Unconjugate

技术规格

Background:	This gene encodes the second largest subunit of RNA polymerase III, the polymerase responsible for synthesizing transfer and small ribosomal RNAs in eukaryotes. The largest subunit and the encoded protein form the catalytic center of RNA polymerase III. Mutations in this gene are a cause of hypomyelinating leukodystrophy. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
Applications:	WB
Name of antibody:	POLR3B
Immunogen:	Synthesized peptide derived from internal of human RPC2.

Full name:	polymerase (RNA) III (DNA directed) polypeptide B
Synonyms :	C128; HLD8; RPC2; INMAP
SwissProt:	Q9NW08
WB Predicted band size:	128 kDa
WB Positive control:	LOVO cells lysate
WB Recommended dilution:	500-3000

