

兔抗 SLC24A1 多克隆抗体

中文名称：兔抗 SLC24A1 多克隆抗体

英文名称：Anti-SLC24A1 rabbit polyclonal antibody

别名：solute carrier family 24 member 1; NCKX; RODX; NCKX1; CSNB1D; HsT17412

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SLC24A1

反应种属：Human, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the potassium-dependent sodium/calcium exchanger protein family. The encoded protein plays an important role in sodium/calcium exchange in retinal rod and cone photoreceptors by mediating the extrusion of one calcium ion and one potassium ion in exchange for four sodium ions. Mutations in this gene may play a role in congenital stationary night blindness. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
Applications:	ELISA, IHC

Name of antibody:	SLC24A1
Immunogen:	Synthetic peptide of human SLC24A1
Full name:	solute carrier family 24 member 1
Synonyms:	NCKX; RODX; NCKX1; CSNB1D; HsT17412
SwissProt:	O60721
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-100

