

兔抗 KLHL4 多克隆抗体

中文名称： 兔抗 KLHL4 多克隆抗体

英文名称： Anti-KLHL4 rabbit polyclonal antibody

别名： KHL4; DKELCHL

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： KLHL4

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the kelch family of proteins, which are characterized by kelch repeat motifs and a POZ/BTB protein-binding domain. It is thought that kelch repeats are actin binding domains. However, the specific function of this protein has not been determined. Alternative splicing of this gene results in two transcript variants encoding different isoforms.
Applications:	ELISA, IHC
Name of antibody:	KLHL4
Immunogen:	Synthetic peptide of human KLHL4
Full name:	kelch like family member 4

Synonyms:	KHL4; DKELCHL
SwissProt:	Q9C0H6
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	25-100

