

兔抗 KLHL8 多克隆抗体

中文名称: 兔抗 KLHL8 多克隆抗体

英文名称: Anti-KLHL8 rabbit polyclonal antibody

相关类别: 一抗

储 存: 冷冻 (-20℃)

宿 主: Rabbit

抗 原: KLHL8

反应种属: Human

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Kelch-like protein 8 is a protein that in humans is encoded by the KLHL8 gene. Substrate-specific adapter of a BCR (BT B-CUL3-RBX1) E3 ubiquitin ligase complex required for The BCR(KLHL8) ubiquitin ligase complex mediates ubiquitination and degradation of RAPSN.
Applications:	ELISA, WB, IHC
Name of antibody:	KLHL8
Immunogen:	Synthetic peptide of human KLHL8
Full name:	kelch-like family member 8
SwissProt:	Q9P2G9
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human colon cancer and Human liver cancer

IHC Recommend dilution:	50-200
WB Predicted band size:	69 kDa
WB Positive control:	Human placenta tissue
WB Recommended dilution:	200-1000



