



## 兔抗 POLR3F 多克隆抗体

中文名称： 兔抗 POLR3F 多克隆抗体

英文名称： Anti-POLR3F rabbit polyclonal antibody

别 名： RNA polymerase III subunit F; RPC6; RPC39

相关类别： 一抗

储 存： 冷冻 (-20°C)

宿 主： Rabbit

抗 原： POLR3F

反应种属： Human, Mouse, Rat


标 记 物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is one of more than a dozen subunits forming eukaryotic RNA polymerase III (RNA Pol III), which transcribes 5S ribosomal RNA and tRNA genes. This protein has been shown to bind both TFIIIB90 and TBP, two subunits of RNA polymerase III transcription initiation factor IIIB (TFIIB). Unlike most of the other RNA Pol III subunits, the encoded protein is unique to this polymerase. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	POLR3F

<b>Immunogen:</b>	Fusion protein of human POLR3F
<b>Full name:</b>	RNA polymerase III subunit F
<b>Synonyms:</b>	RPC6; RPC39
<b>SwissProt:</b>	Q9H1D9
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer and Human brain
<b>IHC Recommend dilution:</b>	50-300
<b>WB Predicted band size:</b>	36 kDa
<b>WB Positive control:</b>	PC-3 and HepG2 cell lysates
<b>WB Recommended dilution:</b>	500-2000



**mlbio** 海联生物  
Good elisakit producers

