

兔抗 SLC25A2 多克隆抗体

中文名称：兔抗 SLC25A2 多克隆抗体

英文名称： Anti-SLC25A2 rabbit polyclonal antibody

别名： ORC2; ORNT2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SLC25A2

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Located between the protocadherin beta and gamma gene clusters on chromosome 5, this intronless gene encodes a protein that is highly similar to an ornithine transporter localized in the mitochondrial inner membrane. The encoded protein most likely plays a role in metabolism as a mitochondrial transport protein.
Applications:	ELISA, IHC
Name of antibody:	SLC25A2
Immunogen:	Fusion protein of human SLC25A2
Full name:	solute carrier family 25 (mitochondrial carrier; ornithine)

	ine transporter) member 2
Synonyms :	ORC2; ORNT2
SwissProt:	Q9BXI2
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	20-100

