

## 兔抗 ZFP36L1(Phospho-Ser92) 多克隆抗体

中文名称：兔抗 ZFP36L1(Phospho-Ser92) 多克隆抗体

英文名称：Anti-ZFP36L1(Phospho-Ser92) rabbit polyclonal antibody

别名：BRF1; ERF1; cMG1; ERF-1; Berg36; TIS11B; RNF162B

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：ZFP36L1(Phospho-Ser92)

反应种属：Human Mouse

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

#### Background:

This gene is a member of the TIS11 family of early response genes, which are induced by various agonists such as the phorbol ester TPA and the polypeptide mitogen EGF. This gene is well conserved across species and has a promoter that contains motifs seen in other early-response genes. The encoded protein contains a distinguishing putative zinc finger domain with a repeating cys-his motif. This putative nuclear transcription factor most likely functions in regulating the response to growth factors. Alternatively spliced transcript variants enc

	oding different isoforms have been found for this gene.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	ZFP36L1(Phospho-Ser92)
<b>Immunogen:</b>	Peptide sequence around phosphorylation site of Serine 92(S-F-S(p)-E-G) derived from Human TISB.
<b>Full name:</b>	ZFP36 ring finger protein-like 1
<b>Synonyms :</b>	BRF1; ERF1; cMG1; ERF-1; Berg36; TIS11B; RNF162 B
<b>SwissProt:</b>	Q07352
<b>WB Predicted band size:</b>	36 kDa
<b>WB Positive control:</b>	HepG2 cells lysates
<b>WB Recommended dilution:</b>	500-1000

