

兔抗 KLK1 多克隆抗体

- 中文名称：兔抗 KLK1 多克隆抗体
- 英文名称：Anti-KLK1 rabbit polyclonal antibody
- 别名：hK1; KLKR; Klk6
- 抗原：KLK1
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 相关类别：一抗
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	Kallikreins are a subgroup of serine proteases having diverse physiological functions. Growing evidence suggests that many kallikreins are implicated in carcinogenesis and some have potential as novel cancer and other disease biomarkers. This gene is one of the fifteen kallikrein subfamily members located in a cluster on chromosome 19. This protein is functionally conserved in its capacity to release the vasoactive peptide, Lys-bradykinin, from low molecular weight kininogen.
Applications:	ELISA, IHC
Name of antibody:	KLK1

Immunogen:	Fusion protein of human KLK1
Full name:	kallikrein 1
Synonyms :	hK1; KLKR; Klk6
SwissProt:	P06870
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human prostate cancer and human liver cancer
IHC Recommend dilution:	50-200

