

兔抗 ZFP36L2 多克隆抗体

中文名称：兔抗 ZFP36L2 多克隆抗体

英文名称： Anti-ZFP36L2 rabbit polyclonal antibody

别名： BRF2; ERF2; ERF-2; TIS11D; RNF162C

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ZFP36L2

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is a member of the TIS11 family of early response genes. Family members are induced by various agonists such as the phorbol ester TPA and the polypeptide mitogen EGF. The encoded protein contains a distinguishing putative zinc finger domain with a repeating cys-his motif. This putative nuclear transcription factor most likely functions in regulating the response to growth factors.
Applications:	ELISA, IHC
Name of antibody:	ZFP36L2
Immunogen:	Synthetic peptide of human ZFP36L2
Full name:	ZFP36 ring finger protein-like 2

Synonyms :	BRF2; ERF2; ERF-2; TIS11D; RNF162C
SwissProt:	P47974
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human esophagus cancer and Human brain
IHC Recommend dilution:	25-100

