

兔抗 SLC25A26 多克隆抗体

中文名称：兔抗 SLC25A26 多克隆抗体

英文名称： Anti-SLC25A26 rabbit ployclonal antibody

别名： SAMC

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SLC25A26

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Mitochondrial carriers, including SLC25A26, are a family of transport proteins found mostly in the inner membranes of mitochondria. They shuttle metabolites and cofactors through the mitochondrial membrane
Applications:	ELISA, IHC
Name of antibody:	SLC25A26
Immunogen:	Fusion protein of human SLC25A26
Full name:	solute carrier family 25 (S-adenosylmethionine carrier), member 26
Synonyms :	SAMC

SwissProt:	Q70HW3
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

