

## 兔抗 KLK3 多克隆抗体

中文名称：兔抗 KLK3 多克隆抗体

英文名称：Anti-KLK3 rabbit polyclonal antibody

别名：APS, PSA

相关类别：一抗

储存：冷冻（-20℃）

抗原：KLK3

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	KLK3 is a member of the kallikrein-related peptidase family secreted by the epithelial cells of the prostate gland and lead to the liquefaction of the seminal coagulum. It is thought to function normally in the liquefaction of seminal coagulum, presumably by hydrolysis of the high molecular mass seminal vesicle protein. Serum level of this protein, called PSA in the clinical setting, is useful in the diagnosis and monitoring of prostatic carcinoma.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	KLK3
<b>Immunogen:</b>	Fusion protein of human KLK3

<b>Full name:</b>	Kallikrein-3
<b>Synonyms :</b>	APS, PSA
<b>SwissProt:</b>	P07288
<b>ELISA Recommended dilution:</b>	500-2000
<b>IHC positive control:</b>	Human cervical cancer and human lung cancer
<b>IHC Recommend dilution:</b>	25-50

