

兔抗 SLC25A38 多克隆抗体

中文名称：兔抗 SLC25A38 多克隆抗体

英文名称： Anti-SLC25A38 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： SLC25A38

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is a member of the mitochondrial carrier family. The encoded protein is required during erythropoiesis and is important for the biosynthesis of heme. Mutations in this gene are the cause of autosomal congenital sideroblastic anemia.
Applications:	ELISA, IHC
Name of antibody:	SLC25A38
Immunogen:	Synthetic peptide of human SLC25A38
Full name:	solute carrier family 25, member 38
SwissProt:	Q96DW6
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human breast cancer
IHC Recommend dilution:	25-100

