

## 兔抗 PONC 多克隆抗体

中文名称：兔抗 PONC 多克隆抗体

英文名称：Anti-PONC rabbit polyclonal antibody

别名：PPNOC

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PONC

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Nociceptin is the ligand of the opioid receptor-like receptor (OPRL1). It may act as a transmitter in the brain by modulating nociceptive and locomotor behavior. May be involved in neuronal differentiation and development
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	PONC
<b>Immunogen:</b>	Synthetic peptide of human PONC
<b>Full name:</b>	prepronociceptin
<b>Synonyms :</b>	PPNOC
<b>SwissProt:</b>	Q13519
<b>ELISA Recommended dilution:</b>	2000-10000

<b>WB Predicted band size:</b>	20 kDa
<b>WB Positive control:</b>	Human liver cancer tissue
<b>WB Recommended dilution:</b>	1000-5000

