

兔抗 POP4 多克隆抗体

- 中文名称：兔抗 POP4 多克隆抗体
- 英文名称：Anti-POP4 rabbit polyclonal antibody
- 别名：POP4 homolog, ribonuclease P/MRP subunit; RPP29
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：POP4
- 反应种属：Human, Mouse, Rat
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes one of the protein subunits of the small nucleolar ribonucleoprotein complexes: the endoribonuclease for mitochondrial RNA processing complex and the ribonuclease P complex. The encoded protein is localized to the nucleus and associates directly with the RNA component of these complexes. This protein is involved in processing of precursor RNAs. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	POP4
Immunogen:	Fusion protein of human POP4

Full name:	POP4 homolog, ribonuclease P/MRP subunit
Synonyms:	RPP29
SwissProt:	O95707
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	50-300

